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http://site.geht.org/UserFiles/file/tableau_synthese_GEHT2015.pdf

Bibliographie

1. Adcock Funk, D.M., et al., *Collection, transport and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays; approved guideline*, in *CLSI document H21-A5*, P.A. Wayne and C.L.S. Institute, Editors. 2008.
2. Mackie, I., et al., *Guidelines on the laboratory aspects of assays used in haemostasis and thrombosis*. *Int J Lab Hematol*, 2012. **35**(1): p. 1-13.
3. Pengo, V., et al., *Update of the guidelines for lupus anticoagulant detection. Subcommittee on Lupus Anticoagulant/Antiphospholipid Antibody of the Scientific and Standardisation Committee of the International Society on Thrombosis and Haemostasis*. *J Thromb Haemost*, 2009. **7**(10): p. 1737-40.
4. Moore, G.W., *Commonalities and contrasts in recent guidelines for lupus anticoagulant detection*. *Int J Lab Hematol*, 2014. **36**(3): p. 364-73.
5. Ledford-Kraemer, M.R., et al., *Laboratory testing for the lupus anticoagulant; Approved guideline - first edition*, in *CLSI document H60-A*, P.A. Wayne and C.L.S. Institute, Editors. 2014.
6. Adcock Funk, D.M., G. Lippi, and E.J. Favaloro, *Quality standards for sample processing, transportation, and storage in hemostasis testing*. *Semin Thromb Hemost*, 2012. **38**(6): p. 576-85.
7. Brien, W.F., et al., *Lupus anticoagulant testing: effect of the platelet count on the activated partial thromboplastin time*. *Br J Biomed Sci*, 1993. **50**(2): p. 114-6.
8. Tripodi, A., et al., *Standardization of activated protein C resistance testing: effect of residual platelets in frozen plasmas assessed by commercial and home-made methods*. *Br J Haematol*, 2003. **120**(5): p. 825-8.
9. Suchsland, J., et al., *Optimizing centrifugation of coagulation samples in laboratory automation*. *Clin Chem Lab Med*, 2014. **52**(8): p. 1187-91.
10. Carroll, W.E., et al., *The significance of platelet counts in coagulation studies*. *J Med*, 2001. **32**(1-2): p. 83-96.
11. Barnes, P.W., C.S. Eby, and M. Lukoszyk, *Residual platelet counts in plasma prepared for routine coagulation testing with the Beckman Coulter Power processor*. *Laboratory Hematology*, 2002. **8**: p. 205-9.
12. Lippi, G., et al., *Influence of the centrifuge time of primary plasma tubes on routine coagulation testing*. *Blood Coagul Fibrinolysis*, 2007. **18**(5): p. 525-8.
13. Nelson, S., A. Pritt, and R.A. Marlar, *Rapid preparation of plasma for 'Stat' coagulation testing*. *Arch Pathol Lab Med*, 1994. **118**(2): p. 175-6.
14. Pappas, A.A., et al., *Rapid preparation of plasma for coagulation testing*. *Arch Pathol Lab Med*, 1991. **115**(8): p. 816-7.
15. Boudaoud, L., et al., *[Rapid centrifugation for routine coagulation testing]*. *Ann Biol Clin (Paris)*, 2006. **64**(4): p. 315-7.
16. Sultan, A., *Five-minute preparation of platelet-poor plasma for routine coagulation testing*. *East Mediterr Health J*, 2010. **16**(2): p. 233-6.
17. Lippi, G., et al., *Influence of residual platelet count on routine coagulation, factor VIII, and factor IX testing in postfreeze-thaw samples*. *Semin Thromb Hemost*, 2013. **39**(7): p. 834-9.
18. van Geest-Daelderop, J.H., et al., *Preanalytical variables and off-site blood collection: influences on the results of the prothrombin time/international normalized ratio test and implications for monitoring of oral anticoagulant therapy*. *Clin Chem*, 2005. **51**(3): p. 561-8.

19. Lippi, G., et al., *Influence of centrifuge temperature on routine coagulation testing*. Clin Chem, 2006. **52**(3): p. 537-8.
20. Cofrac, *Recommandations concernant la surveillance des centrifugeuses*. 2014.
21. Daves, M., et al., *Influence of centrifuge brake on residual platelet count and routine coagulation tests in citrated plasma*. Blood Coagul Fibrinolysis, 2014. **25**(3): p. 292-5.
22. Favaloro, E.J., *Preanalytical variables in coagulation testing*. Blood Coagul Fibrinolysis, 2007. **18**(1): p. 86-9.
23. Lippi, G., et al., *Dishomogeneous separation of citrated plasma in primary collection tubes for routine coagulation testing*. Blood Coagul Fibrinolysis, 2008. **19**(4): p. 330-2.
24. Da Rin, G. and G. Lippi, *Total Laboratory Automation of Routine Hemostasis Testing*. J Lab Autom, 2014. **19**(4): p. 419-422.
25. Flamant, F., et al., *[Comparison between the centrifugation on MPA C10 (Roche Diagnostics) and the centrifugation according recommendations of GEHT (Groupe d'etude de l'hemostase et de la thrombose) for the daily hemostasis assays]*. Ann Biol Clin (Paris), 2014. **72**(2): p. 231-5.
26. Sedille-Mostafaie, N., et al., *Advancing haemostasis automation--successful implementation of robotic centrifugation and sample processing in a tertiary service hospital*. Clin Chem Lab Med, 2013. **51**(6): p. 1273-8.